

Anti-LILRA5 antibody

Catalog Number: 176436

Product name

Polyclonal

Anti-LILRA5 antibody

Ig Type**Specificity**

IgG

Human

Applications**Antibody description**

WB

Rabbit polyclonal antibody to LILRA5

Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human lilra5. aa range:141-190

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality

Fig1.; Western Blot analysis of HeLa cells using CD85f Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000