

Anti-DDIAS antibody



Catalog Number: 176444

Product name

Anti-DDIAS antibody

Ig Type

IgG

Specificity

Human

Applications

WB

Antibody description

Rabbit polyclonal antibody to DDIAS

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human noxin. aa range:111-160

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1.; Western Blot analysis of 823, 293T, AD293, Hela, HepG2 cells using NOXIN Polyclonal Antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000