

Anti-PPP1R15B antibody



Catalog Number: 176445

Product name

Anti-PPP1R15B antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PPP1R15B

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human ppp1r15b. aa range:411-460

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations





Fig1:: Western Blot analysis of NIH-3T3 cells using PPP1R15B Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at



1:20000

Fig2:: Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100