

Anti-CD96 antibody



Catalog Number: 176447

Product name

Anti-CD96 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CD96

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human cd96. aa range:291-340

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 293-UV cells using CD96 Polyclonal Antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000