

# Anti-IGLL1 antibody

Catalog Number: 176451

## Product name

Anti-IGLL1 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to IGLL1

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human igll1. aa range:151-200

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1.; Western Blot analysis of KB cells using Lambda 5 Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at




Fig2.; Immunohistochemical analysis of paraffin-embedded human-lymph, antibody was diluted at 1:100