

# Anti-TNS2 antibody



Catalog Number: 176455

## Product name

Anti-TNS2 antibody

## Ig Type

IgG

## Specificity

Human, Mouse

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to TNS2


## Dilutions

WB: 1:500-1:2000

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tensin-2. aa range:439-488

## Validations



## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

Fig1:; Western blot analysis of A549 3T3 using Tensin-2 antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000