

Anti-ACER1 antibody



Catalog Number: 176459

Product name

Anti-ACER1 antibody

Ig Type

IgG

Specificity

Human

Applications

WB

Antibody description

Rabbit polyclonal antibody to ACER1

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was synthesized peptide derived from asah3 . at aa range: 100-180

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal




Fig1.; Western blot analysis of HepG2 using ASAH3 antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000