

Anti-TNFRSF17 antibody



Catalog Number: 176465

Product name

Anti-TNFRSF17 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to TNFRSF17

Preparation

This antigen of this antibody was synthesized peptide derived from the internal region of human cd269.

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of L929 293T MCF7 HELA using CD269 antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000