

# Anti-TNFRSF17 antibody



Catalog Number: 176465

## General Information

### Product name:

Anti-TNFRSF17 antibody

### Immunogen:

### Applications:

WB

### Specificity:

Human

### Antibody description:

Rabbit polyclonal antibody to TNFRSF17

## Properties

### Preparation:

This antigen of this antibody was synthesized peptide derived from the internal region of human cd269.

### Formulation:

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

### Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

### Clonality:

Polyclonal

### Ig Type: IgG

## Applications

WB: 1:500-1:2000

## Validations




Fig1.; Western blot analysis of L929 293T MCF7 HELA using CD269 antibody. Antibody was diluted at 1:2000. Secondary antibody (catalog# HA1001) was diluted at 1:20000