

# Anti-LAIR2 antibody



Catalog Number: 176466

## Product name

Polyclonal

Anti-LAIR2 antibody

## Ig Type

## Specificity

IgG

Human

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to LAIR2

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human lair2. aa range:21-70

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

## Validations




Fig1.; Western blot analysis of RAT-KIDNEY lysis using LAIR2 antibody. Antibody was diluted at 1:1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000