

# Anti-TP53INP2 antibody



Catalog Number: 176469

## Product name

Anti-TP53INP2 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to TP53INP2

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human tp53inp2. aa range:161-210

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1.; Western blot analysis of SH-SY5Y 293T mouse-kidney lysis using TP53INP2 antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).