

Anti-CLEC4A antibody

Catalog Number: 176473

Product name

Anti-CLEC4A antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to CLEC4A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human clec4a. aa range:161-210

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of HEPG2 cells using CLECSF6 Polyclonal Antibody diluted at 1:500. Secondary antibody (catalog# HA1001) was diluted at 1:20000




Fig2.; Immunohistochemical analysis of paraffin-embedded human-colon-cancer, antibody was diluted at 1:200