

Anti-AVPR1B antibody

Catalog Number: 176494

Product name

Anti-AVPR1B antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to AVPR1B

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human avpr1b. aa range:271-320

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

IHC-P

Dilutions

IHC-P: 1:50-200


Validations

Fig1.; Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200