

Anti-AVPR1B antibody



Catalog Number: 176494

Product name

Polyclonal

Anti-AVPR1B antibody

Ig Type**Specificity**

IgG

Human, Mouse, Rat

Applications**Antibody description**

IHC-P

Rabbit polyclonal antibody to AVPR1B

Dilutions**Preparation**

IHC-P: 1:50-200

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human avpr1b. aa range:271-320

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality**Validations**

Fig1.; Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:200