

Anti-APLF antibody



Catalog Number: 176506

Product name

Anti-APLF antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to APLF

Preparation

This antigen of this antibody was synthesized peptide derived from human aplf polyclonal

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody [catalog#: HA1001) was diluted at 1:20000