

# Anti-CHURC1-FNTB antibody



Catalog Number: 176522

## Product name

Anti-CHURC1-FNTB antibody

## Specificity

Human, Rat

## Antibody description

Rabbit monoclonal antibody to CHURC1-FNTB

## Preparation

This antigen of this antibody was recombinant protein

## Formulation

Liquid, 1\*TBS (pH7.4), 0.05% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB, ICC/IF

## Dilutions

WB: 1:1,000-1:2,000

ICC/IF: 1:50-1:200

## Validations




Fig1: Western blot analysis of FNTB on K562 cells lysates using anti-FNTB antibody at 1/1,000




Fig2: ICC staining FNTB in JAR cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.