

Anti-CD68 antibody



Catalog Number: 176526

General Information

Product name:

Anti-CD68 antibody

Immunogen:

Applications:

WB

Specificity:

Human

Antibody description:

Rabbit monoclonal antibody to CD68

Properties

Preparation:

This antigen of this antibody was synthetic peptide within human cd68 aa 30-70.

Formulation:

Liquid, 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage:

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality:

Monoclonal

Ig Type: IgG

Applications

WB: 1:500-1:1,000

Validations




Fig1.: Western blot analysis of CD68 on SH-SY5Y cell lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.