

# Anti-MYH7B antibody



Catalog Number: 176590

## Product name

Anti-MYH7B antibody

## Specificity

Rat, Human, Mouse

## Antibody description

Rabbit monoclonal antibody to MYH7B

## Preparation

This antigen of this antibody was recombinant protein within c-terminal human myh7b.

## Formulation

Liquid, 1\*TBS (pH7.4), 0.05% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB:1:500-1:1,000

## Validations




Fig1:; Western blot analysis of MYH7B on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody ( 1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.; Positive control;; Lane 1: rat skeletal muscle tissue lysate; Lane 2: rat heart tissue lysate