

# Anti-Ntrk2 antibody



Catalog Number: 176595

## Product name

Anti-Ntrk2 antibody

## Specificity

Mouse, Rat

## Antibody description

Rabbit monoclonal antibody to Ntrk2

## Preparation

This antigen of this antibody was recombinant protein within human trkb aa 100-450.

## Formulation

Liquid, 1\*TBS (pH7.4), 0.05% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2,000

## Validations




Fig1.; Western blot analysis of TrkB on different lysates. Proteins were transferred to a PVDF membrane and blocked with 5% BSA in PBS for 1 hour at room temperature. The primary antibody (1/500) was used in 5% BSA at room temperature for 2 hours. Goat Anti-Rabbit IgG - HRP Secondary Antibody (HA1001) at 1:5,000 dilution was used for 1 hour at room temperature.; Positive control;; Lane 1: mouse brain tissue lysate; Lane 2: mouse cerebellum tissue lysate; Lane 3: rat cerebellum tissue lysate