

# Anti-Tmem39a antibody



Catalog Number: 176706

## Product name

Anti-Tmem39a antibody

## Specificity

Human, Mouse

## Antibody description

Rabbit polyclonal antibody to Tmem39a

## Preparation

This antigen of this antibody was peptide

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

ICC, IHC-P, FC

## Dilutions

ICC: 1:50-1:200

IHC-P: 1:50-1:200

FC: 1:50-1:100

## Validations




Fig1: Immunocytochemical staining of HeLa cells using anti-TMEM39a rabbit polyclonal antibody.





Fig2: Immunocytochemical staining of PANC-1 cells using anti-TMEM39a rabbit polyclonal



antibody.  
Fig3: Immunohistochemical analysis of paraffin-embedded human pancreas tissue using anti-TMEM39a rabbit polyclonal antibody.

# Anti-Tmem39a antibody

Catalog Number: 176706

---




Fig4: Flow cytometric analysis of HeLa cells with TMEM39a antibody at 1/100 dilution (blue) compared with an unlabelled control (cells without incubation with primary antibody; red). Alexa Fluor 488-conjugated Goat anti rabbit IgG was used as the secondary antibody.