

# Anti-Figla antibody



Catalog Number: 176712

## Product name

Anti-Figla antibody

## Ig Type

IgG

## Specificity

Mouse

## Applications

ICC

## Antibody description

Rabbit polyclonal antibody to Figla

## Dilutions

ICC: 1:100

## Preparation

This antigen of this antibody was synthetic peptide within n-terminal human figla.

## Formulation

Liquid, 1\*PBS (pH7.4), 0.2% BSA, 40% Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Validations




Fig1: Immunocytochemical staining of mouse three-cell embryos using anti-Figla rabbit polyclonal antibody (red). (Green: F-actin, Blue: DAPI)