

Anti-PLN antibody



Catalog Number: 176721

Product name

IgG1

Anti-PLN antibody

Applications

Specificity

WB, IP, IF

Human, Mouse, Rat

Dilutions

Antibody description

WB: 1:100-1:1,000

Mouse monoclonal antibody to PLN

IP: 1-2 µg per 100-500 µg of total protein(1 ml of cell lysate)

Preparation

IF: 1:50-1:500


This antigen of this antibody was specific for an epitope mapping between amino acids 7-32 at the n-terminus of phospholamban of human origin.

Formulation

Validations

Liquid, 1*TBS (pH7.4), 1%BSA. Preservative: 0.05% Sodium Azide.

Storage



Store at +4°C.

Clonality

Fig1: Western blot analysis of phospholamban expression in mouse heart (A) and human heart (B) tissue extracts.

Monoclonal

Ig Type