

# Anti-SKA3 antibody



Catalog Number: 176723

## Product name

Anti-SKA3 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Mouse monoclonal antibody to SKA3

## Preparation

This antigen of this antibody was amino acids 154-166 within an internal region of ska3 of human origin

## Formulation

Liquid, 1\*TBS (pH7.4), 1%BSA, 40%Glycerol.  
Preservative: 0.05% Sodium Azide.

## Storage

Store at +4°C

## Clonality

Monoclonal

## Ig Type

IgG1

## Applications

WB, IP, IF

## Dilutions

WB: 1:100-1:1,000

IP: 1-2 µg per 100-500 µg of total protein(1 ml of cell lysate)

IF: 1:50-1:500

## Validations




Fig1: Western blot analysis of Ska3 expression in non-transfected (A) and human Ska3 transfected (B) 293T whole cell lysates.