

# Eribulin Monoclonal Antibody



Catalog Number: 176729

## Product name

Eribulin Monoclonal Antibody

## Specificity

Chemical

## Antibody description

Product has been tested in direct ELISA.

## Preparation

Eribulin conjugated with BSA. Protein G affinity purified

## Formulation

PBS(pH7.4), lyophilized. **Reconstitution: Adding sterile water, prepare a stock solution of 2.0 mg/ml.**

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1,  $\kappa$


## Applications

ELISA

## Dilutions


ELISA: 1:5,000-1:20,000

## Validations



ELISA Binding Assay of Anti-Eribulin Antibody (MM461B) to different proteins.

ELISA analysis of Anti-Eribulin Monoclonal Antibody with Eribulin-BSA, Eribulin-OVA, BSA and OVA. Anti-Eribulin Monoclonal Antibody can bind to Eribulin-BSA or Eribulin-OVA but not carrier proteins.



ELISA Binding Assay of Anti-Eribulin Antibody (MM461B).

Immobilized Eribulin ADC at 2 $\mu$ g/mL (100 $\mu$ L/well) can bind anti-Eribulin mouse monoclonal antibody.