

DXd Monoclonal Antibody



Catalog Number: 176730

Product name

DXd Monoclonal Antibody

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

DXd conjugated with BSA. Protein G affinity purified

Formulation

PBS(pH7.4)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1, κ

Applications


ELISA

Dilutions

ELISA: 1:5,000-1:20,000

Validations


Mouse anti-DXd Antibody ELISA
0.1 μ g/ml



ELISA Binding Assay of Anti-DXd Antibody (MM438I) to different proteins.

ELISA analysis of Anti-DXd Monoclonal Antibody with Dxd-BSA, BSA and Dxd ADC. Anti-DXd Monoclonal Antibody can bind to Dxd-BSA or Dxd ADC but not carrier protein.

anti-DXd mouse mAb MM438I



ELISA Binding Assay of Anti-DXd Antibody (MM438I)

Immobilized Dxd ADC at 2 μ g/mL (100 μ L/well) can bind anti-DXd mouse monoclonal antibody.