

SN-38 Monoclonal Antibody(MM48A)



Catalog Number: 176732

Product name

SN-38 Monoclonal Antibody(MM48A)

Specificity

Chemical

Antibody description

Product has been tested in direct ELISA.

Preparation

SN-38 conjugated with BSA. Protein G affinity purified

Formulation

PBS(pH7.4)

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG1, κ


Applications

ELISA

Dilutions


ELISA: 1:5,000-1:20,000

Validations



ELISA Binding Assay of Anti-SN-38 Antibody (MM48A) to different proteins.

ELISA analysis of Anti-SN38 Monoclonal Antibody with SN-38 -BSA, SN-38 ADC and BSA. Anti-SN-38 Monoclonal Antibody can bind to SN-38-BSA or SN-38 ADC but not carrier protein.



ELISA Binding Assay of Anti-SN-38 Antibody (MM48A)

Immobilized SN-38 ADC at 2 μ g / mL (100 μ L/well) can bind anti-SN-38 mouse monoclonal antibody.