

Anti-TOMM20 Rabbit antibody



Catalog Number: 167407

Product name

Anti-TOMM20 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

TOMM20 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human TOMM20

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/50


ICC/IF: 1/100

FC: 1/50

Validations



Immunocytochemistry of TOMM20(green) in HeLa cells using TOMM20 Rabbit pAb at dilution 1/50, and DAPI(blue)



and DAPI(blue)

Western blot detection of TOMM20 in MCF7, T47D cell lysates using TOMM20 Rabbit pAb(1:1000 diluted). Predicted band size: 16KDa. Observed band size: 16KDa.