

# CDC42EP4 antibody



Catalog Number: 109107

## Product name

CDC42EP4 antibody

## Clonality

Polyclonal

## Immunogen

[Human CDC42EP4 Recombinant protein \(His tag\)](#)

## Ig Type

Rabbit IgG

## Specificity

human, mouse, rat; other species not tested.

## Applications

ELISA, WB

## Antibody description

CDC42EP4 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, A549 cells, mouse lung tissue. Observed molecular weight by Western-blot: 43-45kd


## Dilutions

Recommended Dilution:  
WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of CDC42EP4 recombinant protein (Accession Number: NM\_012121). Purification method: Antigen affinity purified.

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109107(CDC42EP4 antibody) at dilution of 1:1000

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT