

HNRNPH1 antibody



Catalog Number: 111507

Product name

HNRNPH1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

HNRNPH1 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, mouse embryonic stem cells. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 49 kDa

Preparation

This antibody was obtained by immunization of HNRNPH1 recombinant protein (Accession Number: XM_017009416). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111507(HNRNPH1 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:111507(HNRNPH1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog

HNRNPH1 antibody



Catalog Number: 111507

No:111507(HNRNPH1 antibody) at dilution of 1:50 (under 40x lens)