

MBL2 antibody



Catalog Number: 112507

Product name

MBL2 antibody

Specificity

Human; other species not tested.

Antibody description

MBL2 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue. Positive WB detected in human blood tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of MBL2 recombinant protein (Accession Number: NM_001378373). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:112507(MBL2 antibody) at dilution of 1:1000



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:112507(MBL2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:112507(MBL2 Antibody) at dilution of 1:50 (under 40x lens)