

C16orf63 antibody



Catalog Number: 108607

Product name

C16orf63 antibody

Immunogen

[Human C16orf63 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

C16orf63 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 22 kDa

Preparation

This antibody was obtained by immunization of C16orf63 recombinant protein (Accession Number: NM_144600). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:108607(C16orf63 Antibody) at dilution of 1:600