

# DRD5 antibody



Catalog Number: 110017

## Product name

DRD5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DRD5 Rabbit Polyclonal antibody. Positive FC detected in SH-SY5Y cells. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 50-53 kDa

## Preparation

This antibody was obtained by immunization of DRD5 recombinant protein (Accession Number: NM\_000798). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, FC

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110017(DRD5 antibody) at dilution of 1:500



1X10<sup>6</sup> SH-SY5Y cells were stained with 0.2ug DRD5 antibody (Catalog No:110017, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.