

Anti-PHD3 Rabbit antibody



Catalog Number: 167017

Product name

Anti-PHD3 Rabbit antibody

WB: 1/1000

IHC: 1/50

Specificity

Human, Mouse, Rat

ICC/IF: 1/50

IP: 1/20

Antibody description

PHD3 Rabbit polyclonal antibody

Preparation

Antigen: Recombinant protein of human PHD3

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


WB, IHC-P, ICC/IF, IP

Dilutions

Validations



Western blot detection of PHD3 in Rat lung and Mouse lung lysates using PHD3 rabbit pAb (1:1000 diluted). Predicted band size: 27kDa. Observed band size: 27kDa.



Immunohistochemistry of PHD3 in paraffin-embedded Human colon cancer tissue using PHD3 Rabbit pAb at dilution 1/50