

# PHIP antibody



Catalog Number: 113817

## Product name

PHIP antibody

## Immunogen

[Human PHIP Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

PHIP Rabbit Polyclonal antibody. Positive IP detected in A375 cells. Positive WB detected in A375 cells, HeLa cells. Observed molecular weight by Western-blot: 190-207 kDa

## Preparation

This antibody was obtained by immunization of PHIP recombinant protein (Accession Number: BC008909). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:1000-1:10000

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:113817(PHIP antibody) at dilution of 1:3000



IP Result of anti-PHIP (IP:Catalog No:113817, 4ug; Detection:Catalog No:113817 1:3000) with A375 cells lysate 3600ug.