

Anti-Citrate synthetase Rabbit antibody



Catalog Number: 166817

Product name

Anti-Citrate synthetase Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Citrate synthetase Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Citrate synthetase

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC

Dilutions


WB: 1/1000

IHC: 1/50


ICC/IF: 1/20-1/100

FC: 1/20


Validations



Western blot detection of Citrate synthetase in K562, Rat Brain, C6, 3T3, Hela cell lysates using Citrate synthetase Rabbit pAb(1:1000 diluted). Predicted band size: 52kDa. Observed band size: 45kDa.



Immunohistochemistry of Citrate synthetase in paraffin-embedded Human tonsil using Citrate synthetase Rabbit pAb at dilution 1/50



Immunocytochemistry of Citrate synthetase(green) in Hela cells using Citrate synthetase Rabbit pAb at dilution 1/50, and DAPI(blue)