

# DUSP15 antibody



Catalog Number: 110117

## Product name

DUSP15 antibody

## Immunogen

[Human DUSP15 Recombinant protein \(GST tag & His tag\)](#)

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

DUSP15 Rabbit Polyclonal antibody. Positive WB detected in A375 cells. Observed molecular weight by Western-blot: 32kd

## Preparation

This antibody was obtained by immunization of DUSP15 recombinant protein (Accession Number: NM\_080611). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:1000

## Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110117(DUSP15 antibody) at dilution of 1:300