

AXIN1 antibody



Catalog Number: 108517

Product name

AXIN1 antibody

Immunogen

[Human AXIN1 Recombinant protein \(GST tag & His tag\)](#)

Specificity

Human; other species not tested.

Antibody description

AXIN1 Rabbit Polyclonal antibody. Positive WB detected in HT-1080 cells, A431 cells, A549 cells. Observed molecular weight by Western-blot: 100-110kd

Preparation

This antibody was obtained by immunization of AXIN1 recombinant protein (Accession Number: NM_181050). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



HT-1080 cells were subjected to SDS PAGE followed by western blot with Catalog No:117078(AXIN1 antibody) at dilution of 1:500