

G-CSFR antibody



Catalog Number: 110917

Product name

G-CSFR antibody

Specificity

Human; other species not tested.

Antibody description

G-CSFR Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue. Observed molecular weight by Western-blot: 85kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_000760). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110917(G-CSF-R antibody) at dilution of 1:300