

Anti-Phospho-HSF1 (Ser326) Rabbit antibody



Catalog Number: 167217

Product name

Anti-Phospho-HSF1 (Ser326) Rabbit antibody

FC: 1/200

IP: 1/20

Specificity

Human

Antibody description

Phospho-HSF1 (Ser326) Rabbit polyclonal antibody

Preparation

Antigen: A synthetic phosphopeptide corresponding to residues surrounding Ser326 of human HSF1

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P, ICC/IF, FC, IP


Dilutions

WB: 1/1000-1/5000


IHC: 1/20

ICC/IF: 1/50

Validations




Immunocytochemistry of HSF1 (Phospho-Ser326)(green) in Hela cells using HSF1 (Phospho-Ser326) Rabbit pAb at dilution 1/50, and



DAPI(blue)

Immunohistochemistry of HSF1 (Phospho-Ser326) in paraffin-embedded Human breast cancer tissue using HSF1 (Phospho-Ser326) Rabbit pAb at



dilution 1/20

Western blot detection of Phospho-HSF1 (Ser326) in K562, Hela cell lysates using Phospho-HSF1 (Ser326) Rabbit pAb(1:1000 diluted). Predicted band size: 57kDa. Observed band size: 82kDa.