

# TBP antibody



Catalog Number: 107617

## Product name

TBP antibody

## Immunogen

[Human TBP Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat, Pig; other species not tested.

## Antibody description

TBP Mouse Monoclonal antibody. Positive IHC detected in human breast cancer tissue, human liver cancer tissue. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, pig liver tissue. Observed molecular weight by Western-blot: 40 kDa

## Preparation

This antibody was obtained by immunization of TBP recombinant protein (Accession Number: NM\_003194). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IP, IHC

## Dilutions


Recommended Dilution:

WB: 1:2000-1:20000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



Western blot analysis of TBP in various cell lines using Proteintech antibody Catalog No:117338 at a dilution of 1:10000.




IP Result of anti-TBP (IP:Catalog No:117338, 3ug; Detection:Catalog No:117338 1:800) with HEK-293 cells lysate 1800ug.

# TBP antibody

Catalog Number: 107617

---



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:117338(TBP Antibody) at dilution of 1:400 (under 40x lens).