

HuA/B/C/D antibody



Catalog Number: 111717

Product name

HuA/B/C/D antibody

WB: 1:1000-1:10000

IP: 1:200-1:1000

Specificity


Human, Mouse, Rat, Zebrafish; other species not tested.

IHC: 1:20-1:200

Antibody description

HuA/B/C/D Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in human cerebellum tissue, mouse cerebellum tissue, PC-13 cells. Positive IHC detected in human breast cancer tissue. Observed molecular weight by Western-blot: 40 kDa

Validations




Preparation

This antibody was obtained by immunization of HuA/B/C/D recombinant protein (Accession Number: NM_001144774). Purification method: Antigen affinity purified.

human cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:111717(ELAVL4 antibody) at dilution of 1:800

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:111717(ELAVL4 antibody) at dilution of 1:100 (under 10x lens)

Ig Type


Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:



HuA/B/C/D antibody



Catalog Number: 111717

IP Result of anti-ELAVL4 (IP:Catalog No:111717, 4ug; Detection:Catalog No:111717 1:300) with mouse brain tissue lysate 2560ug.