

Cysteine dioxygenase type 1 antibody



Catalog Number: 109717

Product name

Cysteine dioxygenase type 1 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200

Antibody description

Cysteine dioxygenase type 1 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, rat liver tissue. Positive IP detected in mouse liver tissue. Positive IHC detected in human pancreas cancer tissue, human liver cancer tissue, human placenta tissue. Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of Cysteine dioxygenase type 1 recombinant protein (Accession Number: NM_001801). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP


Dilutions

Recommended Dilution:


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109717(CDO1 antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:109717(CDO1 antibody) at dilution of 1:100 (under 10x lens)



Cysteine dioxygenase type 1 antibody



Catalog Number: 109717

IP Result of anti-CDO1 (IP:Catalog No:109717,

3ug; Detection:Catalog No:109717 1:600) with mouse liver tissue lysate 4000ug.