

# Anti-CD229 Rabbit antibody



Catalog Number: 168027

## Product name

Rabbit IgG

Anti-CD229 Rabbit antibody

## Applications

## Specificity

WB, IHC-P

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/1000

CD229 Rabbit monoclonal antibody

IHC: 1/50

## Preparation

Antigen: Recombinant protein of human CD229

## Validations

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA


## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type



Western blot detection of CD229 in HeLa,A549,HL-60,U251 cell lysates using CD229 Rabbit mAb(1:1000 diluted).Predicted band size:72kDa.Observed band size:72kDa.