

APOBEC3G antibody



Catalog Number: 108027

Product name

APOBEC3G antibody

IHC: 1:20-1:200

IF: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

APOBEC3G Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue, human colon cancer tissue. Positive IF detected in MCF-7 cells. Positive WB detected in mouse colon tissue. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of APOBEC3G recombinant protein (Accession Number: NM_021822). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:108027(APOBEC3G Antibody) at dilution of



1:100


Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108027(APOBEC3G Antibody) at dilution of 1:200 (under 10x lens).

Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:108027(APOBEC3G Antibody) at dilution of 1:200 (under 40x lens).

APOBEC3G antibody



Catalog Number: 108027



Immunofluorescent analysis of MCF-7 cells using Catalog No:108027(APOBEC3G Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)