

Anti-Ferritin Heavy Chain Rabbit antibody



Catalog Number: 167427

Product name

Anti-Ferritin Heavy Chain Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Ferritin Heavy Chain Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Ferritin

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions


WB: 1/1000

ICC/IF: 1/20

Validations




Western blot detection of Ferritin in Raw264.7 cell lysates using Ferritin Rabbit pAb(1:1000 diluted). Predicted band size: 21 kDa. Observed



band size: 21 kDa.

Western blot detection of Ferritin in K562, CHO-K1, HeLa cell lysates using Ferritin Rabbit pAb(1:1000 diluted). Predicted band size: 21 kDa. Observed band size: 21 kDa.



Immunocytochemistry of Ferritin (green) in HeLa cells using Ferritin Rabbit pAb at dilution 1/50, and DAPI (blue)