

MRPL38 antibody



Catalog Number: 112827

Product name

MRPL38 antibody

IHC: 1:20-1:200

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

MRPL38 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, human placenta tissue, mouse liver tissue. Positive IF detected in HepG2 cells. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 41kd

Preparation

This antibody was obtained by immunization of MRPL38 recombinant protein (Accession Number: BC013311). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112827(MRPL38 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:112827(MRPL38 antibody) at dilution of 1:100 (under 10x lens)




Immunohistochemical of paraffin-embedded human liver using Catalog No:112827(MRPL38

MRPL38 antibody



Catalog Number: 112827

antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using MRPL38 antibody Catalog No:112827 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).