

Tim23 antibody



Catalog Number: 116127

Product name

Tim23 antibody

Recommended Dilution:

WB: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.


IP: 1:500-1:5000

IHC: 1:20-1:200

Antibody description

Tim23 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue, human heart tissue, human testis tissue. Positive IP detected in mouse heart tissue. Positive WB detected in Jurkat cells, COLO 320 cells, mouse heart tissue, mouse skeletal muscle tissue, rat heart tissue. Observed molecular weight by Western-blot: 22 kDa

Validations



Preparation

This antibody was obtained by immunization of Tim23 recombinant protein (Accession Number: BC062707). Purification method: Antigen affinity purified.

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116127(TIMM23 antibody) at dilution of 1:600

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP


Dilutions



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:116127(TIMM23 antibody) at dilution of 1:50 (under 10x lens)

Tim23 antibody

Catalog Number: 116127



IP Result of anti-Tim23 (IP:Catalog No:116127, 4ug; Detection:Catalog No:116127 1:1000) with mouse heart tissue lysate 3200ug.