

Anti-NM23A (1H9) Mouse antibody



Catalog Number: 168127

Product name

Anti-NM23A (1H9) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

NM23A (1H9) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Synthetic peptide corresponding to human NME1/NDKA protein

Formulation

PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of NME1/NDKA in A549, MCF7, 3T3, C6 and HeLa cell lysates using NME1/NDKA mouse mAb(dilution 1:1000). Predicted band size: 16, 18 kDa. Observed band size: 16, 18 kDa.