

AP4M1 antibody



Catalog Number: 108127

Product name

AP4M1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

AP4M1 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue. Positive WB detected in mouse heart tissue, A431 cells, HeLa cells, human heart tissue, human lung tissue, mouse lung tissue. Observed molecular weight by Western-blot: 50-55 kDa

Preparation

This antibody was obtained by immunization of AP4M1 recombinant protein (Accession Number: XM_047421026). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:108127(AP4M1 Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:108127(AP4M1 antibody) at dilution of 1:50 (under 10x lens)