

EHHADH antibody



Catalog Number: 110227

Product name

EHHADH antibody

Clonality

Polyclonal

Immunogen

[Human EHHADH Recombinant protein \(His tag\)](#)

Ig Type

Rabbit IgG

Specificity

Human; other species not tested.

Applications

ELISA, WB

Antibody description

EHHADH Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue, SMMC-7721 cells. Observed molecular weight by Western-blot: 69 kDa

Dilutions

Recommended Dilution:
WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of EHHADH recombinant protein (Accession Number: NM_001966). Purification method: Antigen affinity purified.


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110227(EHHADH Antibody) at dilution of 1:2000